



Wheelchair



Compact Transport Wheelchair

The Compact Transport Aluminium Wheelchair from Aidapt features an extremely compact and lightweight aluminium frame. Ideal for the occasional or first time user who wants a strong, reliable and easy to use wheelchair that will go anywhere and fold in seconds to go in the boot of your car, coach or aircraft.

- * Ultra light frame
- * Folding backrest (for storage or transport)
- * Fixed armrest
- * Swing away removable leg rests
- * Solid 12" rear castor
- * Care assist hand brake
- * Rear wheel locking brake
- * Weight capacity 115kg



Specifications

Overall height	93	Seat height.....	45
Overall depth	90	Backrest height	42
Overall width	50	Armrest height.....	25
Seat depth.....	40	Unit weight	12
Seat width.....	43	Weight capacity	115

Measurements in cm, weight

Product Code: VA170H Silver
Product Code: VA70B Blue
Product Code: VA170P Pink
Product Code: VA170R Red

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au