



Bedroom

Over Bed Pole

Free standing Over Bed Pole assists in pulling oneself up in bed. Ideally placed under head of the bed.

- * Strong steel frame construction
- * Triangle handle with adjustable strap
- * Rubber suction caps to prevent from slipping
- * Weight Capacity 140kg



Specifications

Overall height	183
Overall depth	80
Overall width	81
Unit weight	11
Weight capacity	140

Measurements in cm, weight in kg

Product Code: RM475

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au