



Fact Sheet

Bathroom *Bariatric*

Over Toilet Aid

The Bariatric 406 Over Toiler Aid constructed from aluminium oval tubing offers extra strength and security in its frame and support. The high density wider pvc seat offers ultimate comfort and comes with easily fitting modesty cover and splash guard. Fitted with large rubber suction tips to reduce slipping.

- * Lightweight toughened powder coated aluminium frame
- * Large suction tips to reduce slipping
- * Wider seat for increased comfort
- * 65cm between arms
- * Modesty cover
- * Splash guard
- * Height adjustable
- * Weight capacity 305kg



Specifications

Overall height 67-79
 Overall depth 54
 Overall width 87
 Seat depth 45
 Seat width 65

Seat height 43-55
 Armrest height 24
 Unit weight 8.7
 Weight capacity 305

Measurements in cm, weight in kg

Product Code: RM406-65

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au