



Wheelchair

Compact Lite Wheelchair

The Compact Lite is as the name suggests compact all round great wheelchair. Foldable, light in weight makes it a great chair not just for the occupant, but also the carer. A smart, sporty and classy looking compact unit, it is comfortable and perfect for every day use.

- * Lightweight frame, easy to manage
- * Swing up leg rests with 2 height options
- * AP - Solid 6" front castor, 12" rear wheel
- * SP - Solid 6" front castor, 20" rear wheel
- * Optional anti-tip available (RM836AT)
- * Rear wheel locking brake
- * Folding backrest
- * Attendant Brake
- * Seat belt included
- * Weight capacity 110kg



Specifications AP/SP

Overall height	92	Seat height.....	46
Overall depth	97/99	Backrest height	44
Overall width	60	Armrest height.....	21
Seat depth.....	38	Unit weight	11/12
Seat width.....	45	Weight capacity	110

Measurements in cm, weight

Product Code: RM836-AP A/P
Product Code: RM836-SP S/P

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au