

Fact Sheet



Wheelchair



Deluxe 20" Self Propelled

The new deluxe lightweight 50cm Self Propelled Aluminium Wheelchair comes with all the features of the VA165 Deluxe lightweight, compact wheelchair, only with extra width in the seat and adjustable armrest height. High grade aluminium with a double cross bar for strength, the lightweight wheelchair also offers ease of movement for the user and quick-folding for trouble-free transportation, with the extra allowance in weight capacity 136kg.

- * Lightweight Self Propelled
- * Care assist hand brake
- * Quick release wheels (for transport)
- * Flip up armrest for side transfer
- * Swing away and removable leg rests
- * Solid tyres (no maintenance required)
- * Can be fitted with elevating leg rests
- * Weight capacity 136kg



Product Code: Silver VAB165S-50

Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Backrest Height	Footrest Height	Unit Weight	Weight Capacity
91	104	69	43	50	53	22-28	39	17	16	136

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au