



GELCOMFORT CUSHIONS



FEATURES

Featuring patented dual 3D Honeycomb and Air-Cell Beta Gel technology, the GelComfort cushion provides a new level of comfort. The dual honeycomb design creates excellent pressure distribution, dispersing pressure on areas prone to discomfort and pressure sores, such as the hips and thighs, and assisting good blood flow.

- Patented technology improves on memory foam and latex thanks to greater stability provided through a cell wall structure, preventing buckling and compression
- Air-Cell Beta Gel keeps the cushion cool and ventilated. Each Air-Cell discharges warm air and

absorbs fresh air with body's movement to ventilate the cushion.

- Beta Gel is made from skin-safe oils, special polymers and gel elastomers, which provide excellent pressure distribution and superior elasticity. Beta Gel protects skin and encourages optimal blood flow through its dispersion of pressure.
- The GelComfort cushion is designed so it can be easily cleaned for continued use. The gel insert can be washed at temperatures up to 50 degrees, while the fabric cover can be machine washed
- Water resistant cover in PU stretch available for Large and Extra Large sizes.

TECHNICAL SPECIFICATIONS

Dimensions

L1: 460 x 420 x 25mm
 L2: 460 x 420 x 50mm
 XL2: 550 x 500 x 50mm

Colours

Grey (L1, L2, XL2), Green (L2), Orange (L2)

Material:

Cover: Polyester 100%
 Gel: Polysiloxanes 100%

Pressure Care Rating:

Medium to High

Order Code

Large (single layer): CM-GC-L1
 Large (dual layer): CM-GC-L2
 Extra Large (dual) : CM-GC-XL2

Warranty

12 months