



Bathroom

Over Toilet Aid - Tool-less Assembly

The RM Compact Over Toilet Aid is a portable easy to assemble shower chair providing secure safe option for those who require a seat in the shower. A very simple click together frame, fixed armrests, clip on and off seat & lid for easy cleaning. Ideal for permanent or temporary solution in the bathroom.



- * Lightweight rust proof aluminium frame
- * Durable powder coat finish
- * Very easy to assemble
- * Dismantle for storage / transportation
- * Height adjustable
- * Complete with splash guard
- * Non-slip rubber tips
- * Weight capacity 110kg



Unassembled
Compact packaging
 58x15x51cm 4kg

Specifications

Overall height	51-74	Seat height.....	40-53
Overall depth	60	Armrest height	19
Overall width	52	Unit weight	2.2
Seat depth	39	Weight capacity.....	110
Seat width.....	47		

Measurements in cm, weight in kg

Product Code: RM390A-KD

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au