



Bathroom *Bariatric*

Shower Stool optional Backrest

The Bariatric 408 Shower Stool constructed from aluminium oval tubing offers extra strength and security in its frame with the option of a backrest for further safety. The PVC seat and backrest provide the ultimate in comfort. Fitted with large rubber suction tips to reduce slipping.

- * Lightweight toughened powder coated aluminium frame
- * Large suction tips to reduce slipping
- * PVC padded seat
- * 65cm between arms
- * Height adjustable
- * Optional PVC padded backrest
- * Weight capacity 305kg



Specifications

Overall height 67-79
 Overall depth 54
 Overall width 87
 Seat depth 35
 Seat width 65

Seat height 43-55
 Armrest height 27
 Unit weight 5.7
 Weight capacity 305

Measurements in cm, weight in kg

Overall height with backrest 86-98
 Backrest width 49

Product Code: RM408-65BR Backrest
Product Code: RM408-65 Stool

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au