



Fact Sheet

Bathroom

Shower Stool - Tool-less Assembly

The RM Compact Aluminium Shower Stool is a portable easy to assemble shower stool providing secure safe option for those who require a seat, without backrest, in the shower. A very simple click together frame; seat already attached, fixed armrests, ideal for permanent or temporary solution in the bathroom.



- * Lightweight rust proof aluminium frame
- * Durable powder coat finish
- * Very easy to assemble
- * Dismantle for storage / transportation
- * Height adjustable
- * Moulded plastic seat with drain holes
- * Non-slip rubber tips
- * Weight capacity 110kg



Unassembled
Compact packaging
65x14x45cm 3kg

Specifications

Overall height..... 66-78.5
Overall depth47
Overall width.....54
Seat depth.....22
Seat width47

Seat height 47-59.5
Armrest height19
Unit weight.....2.2
Weight capacity110

Measurements in cm, weight in kg

Product Code: RM392A-KD

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au