

Bathroom

Aluminium Shower Stool RM403

The RM lightweight aluminium shower stool with durable "S" bend side frame offering maximum comfort and space when seated. With large rubber suction ferrule to prevent slipping and provide maximum safety.

- * Aluminium frame
- * Moulded plastic seat with drainage holes
- * "S" bend side frame provides maximum space
- * Large rubber suction ferrule to prevent slipping
- * Weight capacity 160kg
- * Optional backrest



Product Code: RM403



Product Code: RM403BR

Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Armrest Height	Unit Weight	Weight Capacity
67-82	44	58	33	50	49-64	21	2.8	160

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au